DEFENDER™

Light Curing Opaque Pit & Fissure Sealant

Directions For Use • For Professional Use Only

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DESCRIPTION

Henry Schein Defender is a light cured pit and fissure sealant. It is a very low viscosity resin system with a white opaque color to help identify its retention at recall appointments. It is easily placed and cures quickly to a hard and retentive sealant.

INDICATIONS FOR USE

Henry Schein Defender is used for covering and sealing of pits and fissures in the occlusal surfaces of posterior teeth to aid in the resistance to caries in those surfaces.

FORMULATION

Defender is a dilute BisGMA based resin system with a white opaquer, a microfine silica particle and a photosensitive catalyst system.

DIRECTIONS FOR USE:

1. Clean the teeth with a non-fluoride, oil-free prophyl paste. Wash with an oil-free air/water spray.
2. Isolate teeth with cotton rolls or rubber dam.
3. Blow dry completely with an oil/moisture-free air source.
4. Dispense Etchant (37% phosphoric acid) into a mixing well and apply the etchant to the occlusal grooves of the teeth with a brush. Do not rub against the enamel, but the etchant may be agitated with the brush to increase flow into the fissures. Allow etchant to remain for 30 seconds.
5. Wash area completely for 20 seconds with adequate suction to remove etchant and water. Do not allow saliva to contaminate etched enamel surfaces.
6. Air dry thoroughly using oil/moisture free air.
7. Etched surfaces of teeth should have a dull, white, frosted appearance when etched and dried. If the enamel is not chalky in appearance, repeat the process allowing etchant to remain on the teeth for an additional 30 second increments of time and wash and dry. Several repetitions of this etching procedure may be needed especially on primary teeth and teeth high in fluoride content.
8. Place two drops of Defender into a mixing well. This quantity is sufficient for sealing one quadrant.
9. With a brush apply Defender sealant to each etched, dry occlusal surface taking care to avoid air entrapment. Agitation of the brush bristles may help flow.
10. Hold the end of the light guide as close to the tooth as possible without touching the resin. Expose each coated area to the light curing source for 20-30 seconds. (WARNING: Use a blue curing light (470-480 nanometer with a minimum of 300 mW/cm² intensity).
11. After setting wipe the occlusal surfaces with a cotton roll to remove the oxygen inhibited layer at the surface of the resin, or rinse with copious amounts of water with adequate suction present.
12. Check for complete coverage of the fissures by using an explorer. Incomplete coverage requires re-etching and re-application of sealant. Adjust occlusal surface to avoid occlusal interferences.

WARNINGS

Etchant gels are caustic and can cause tissue damage in contact with skin or eyes. Dentists and assistants should wear gloves and protective eye wear. Patient should also wear eye protection.

If accidental contact occurs with skin, flush immediately with copious amounts of water. In case contact with eyes occur, flush immediately with copious amounts of water and consult a physician.

Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.

Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.

Do not allow saliva or water contamination of unset material.

If a fluoride treatment is planned at the same appointment, apply the fluoride after bonding the sealant, and the sealant is completely set.

STORAGE CONDITIONS

This material is designed to be stored at 35° - 50°F (2° - 10°C). Allow the material to equilibrate to room temperature prior to use. Shelf life is 2 years. Do not store in the proximity of eugenol-containing products. Do not store the material in intense light, or under wet conditions.

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Henry Schein U.K. Holdings Ltd.